Silverbolt™ Instrument Kit

Instruction for Use

STERILE

Rx only

Implants

CE 0086

Choice Spine
400 Erin Drive
Knoxville, TN 37919
USA

Indications for Use:

The Silverbolt™ Screw System Instrument Kit consists of manual surgical instruments intended to be used for implantation of Choice Spine™ orthopedic implants. Prior to use of these instruments, please refer to the product instructions for-use specific to the Silverbolt Screw System and the surgical procedure to be followed.

The Silverbolt Screw System Kit Instruments consists of Guidewires, Target needles, and Tissue Splitter Blades for insertion to assist in dilating the incision and provide access to the bone, and operative site. The Silverbolt Screw System Instrument Kit components are made from materials commonly used in orthopedic procedures.

Precautions:

1. Choice Spine instruments that are labeled, or designed as single use devices should not be re-sterilized, or re-used, and should be disposed of properly after use.
2. Choice Spine instruments should only be used with Choice Spine implants and instruments. Do not attempt to use with other competitive devices.
3. Avoid application of excessive stress on surgical instrumentation.
4. Carefully read and follow any package insert which accompanies the implants to be used with this instrumentation.
5. Carefully inspect all instruments prior to use. Do not use an instrument that is damaged.

Packaging:

The Silverbolt Screw System Instrument Kit is supplied STERILE and intended for single patient use only. DO NOT RESTERILIZE THIS PRODUCT. The sterility can only be assured if the packaging is intact. Do not use this device if the STERILE packaging has been opened or damaged. Contact your local sales representative or distributor for replacement. Remove all packaging material prior to use. Only sterile instruments should be used in surgery.

Patient Education:

It is essential to provide preoperative instructions to the patient. S/he should be made aware of the potential risks of the surgery and the implant limitations. The patient should be instructed to limit postoperative activity, as this will reduce the risk of bent, broken or loose implant components. The patient must be made aware that implant components may bend, break, or loosen even though restrictions in activity are followed.

Storage and Handling:

Implants should be stored in the implant sterilization case in clean, dry, well-ventilated conditions away from floors, ceilings, and outside walls. Store and transport sterile implants in such a way as to maintain sterility and functional integrity. Do not use implants if the sterilization wrap is opened, damaged or wet. Implants should remain covered until needed to avoid contamination. Only those to be implanted should be handled.

Limitations and Restrictions:

Repeated sterilization according to these instructions has a minimal effect on Choice Spine devices. Sterilization equipment varies in performance characteristics and must be validated accordingly. The sterilizing facility is responsible for the routine validation and monitoring of all equipment, materials and personnel used in their facility to ensure the desired results are achieved.

These instructions have been validated as being capable of sterilizing these Choice Spine implants. Any deviations from these procedures must be evaluated for efficacy by the sterilizing facility.

Device Retrieval Efforts:

Should it become necessary to remove any or all of Silverbolt Screw System Instrument Kit components, please call Choice Spine at the number below to receive instructions regarding data collection, including histopathological, mechanical, and adverse event information.

Caution:

Federal Law (USA) restricts this device to sale by or on the order of a physician.

Information:

See www.choicespine.com/patents.html for details.

For product complaints please contact:

Choice Spine, LP
Quality/Regulatory Department
400 Erin Drive
Knoxville, TN 37919
Phone: 865-246-3333; Fax: 865-588-4045

For additional Product information please contact:

Choice Spine, LP
Customer Service Department
400 Erin Drive
Knoxville, TN 37919
Phone: 865-246-3333 or fax: 865-588-4045
customerservice@choicespine.com



Product Complaints:

Any dissatisfaction with the product quality, labeling, or performance should be reported to Choice Spine immediately by the customer or health care provider. Furthermore, Choice Spine should be notified immediately of an implant malfunction by telephone, fax, or written correspondence. When filing a complaint, the name, part number, and lot number of the part should be provided along with the name and address of the person filing the complaint.

Some components may not be currently available. Please contact your Choice Spine representative for additional information. The products discussed herein may be available under different trademarks in different countries. All copyrights, and pending and registered trademarks are property of Choice Spine. For more information on a specific product or trademark, please contact your local Choice Spine representative.

EC	REP
----	-----

MedPass SAS
95 bis Boulevard Pereire,
75017 Paris France
MedPass.ar@medpass.org
Tel: +33 (0) 1 42 12 28 84
Fax: +33 (0) 1 42 12 28 83

Symbol	Definition
	Do not reuse
	Caution, consult instructions for use for warnings and precautions
	Consult instructions for use
	Do not use if package is damaged
	Lot number
	Reference number
	Serial Number
	Sterilized by irradiation
	Use by
	Manufacturer
	Date of Manufacture
	Federal law (USA) restricts this device to sale by or on the order of a physician
	Non-Sterile
	European Medical Devices
	Authorized representative in the European Community

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